

# MULTICUT SOLO

THE MOST EFFICIENT MORCELLATOR FOR YOUR  
UROLOGY & BPH ENUCLEATION TREATMENTS



- ⊕ REDUCED IRRIGATION VOLUME REQUIRED
- ⊕ OPTIMIZED BLADE CUTTING GEOMETRY
- ⊕ MULTIPLE USE BLADES
- ⊕ INCREASED SAFETY

**JENA SURGICAL**  
LASER AT YOUR SIDE



# Technical Specifications

## MultiCut Solo

<b>Device Model</b>	Table-top device
<b>Control Panel</b>	7" color touch screen
<b>Cooling</b>	Active Air Cooling
<b>Accessories</b>	MultiCut Handpiece with two different blades (350 mm and 400 mm), aspiration tube set, waste container
<b>Electrical Requirements</b>	~100-240 V, 50/60 Hz, max. 65-90 VA
<b>Weight</b>	<10 kg
<b>Dimensions</b>	42 cm x 35 cm x 15 cm (W x D x H)



## MultiCut Handpiece

<b>Rotation Speed</b>	500-3,000 rpm
<b>Oscillation Frequency</b>	0 or 2-5 Hz
<b>Morcellator Blades</b>	Blade 350 mm – 15 Fr x 350 mm (diameter x free length), 65 g Blade 400 mm – 15 Fr x 400 mm (diameter x free length), 75 g
<b>Cable Length</b>	3 m
<b>Sterilisation Method</b>	Steam sterilisation
<b>Weight</b>	330 g (without cable and blade)
<b>Dimensions</b>	30 x 200 mm (diameter x free length)

Asclepion Laser Technologies GmbH  
Brüsseler Str. 10 | 07747 Jena | Germany  
[www.jenasurgical.com](http://www.jenasurgical.com)

JenaSurgical is the brand of the surgical business unit of Asclepion Laser Technologies. This brochure is not intended for the market of USA.

